Montage- und Gebrauchsanweisung	5	DEU
Installation and operating instructions	14	ENG

D^{med}® **visano 20**

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1 SAFETY INSTRUCTIONS

1.1 Intended use

The luminaire D^{med}® visano is an examination luminaire. Its intended use is to illuminate the body of the patient to support diagnosis or treatment. An interruption of the diagnosis or treatment caused by a light failure is always possible without a risk for the patient. The luminaire can not be used in operating rooms.

1.2 User profiles

Medical professional

A medical professional is any person who has completed medical training and is working in his or her professional field.

Cleaning professional

A cleaning professional is trained in national and job-related hygiene regulations.

Electrician

An electrician is trained in the fields of electronics and electrical engineering and is familiar with the relevant standards and regulations.

Qualified professional

A qualified professional is capable of mounting and dismounting the luminaire thanks to professional training, knowledge and experience and knowledge of the regulations.

1.3 Safety instructions

- ▶ Operation by a medical professional
- ▶ The instructions form part of the product. They must be stored and made accessible to all subsequent users.
- ▶ Any work on the luminaire (including repairs) must be carried out by qualified electricians only. Mounting must be performed by a qualified professional only.
- ▶ The luminaire must not be altered or manipulated in any way. Only approved original parts must be used. Any use other than the intended use with original parts may give rise to other technical values and life-threatening danger.
- ▶ Operation in potentially explosive areas is prohibited. The luminaire power supply is a potential ignition source.
- ▶ The luminaire must only be operated in dry, dust-free rooms.
- ▶ The luminaire must not be left switched on without supervision.
- ▶ For luminaires in protection class I, it is essential for the earthing conductor to be connected to the luminaire housing.
- ▶ Do not use any luminaire that is damaged. Defective cables are also potential hazards. Do not place the cable close to any heat source or on sharp edges.
- ▶ Damage to the eyes. Never look directly into the light source.
- ▶ Replace any glass that is damaged before operating the luminaire again.
- ▶ In order to prevent electric shocks, only connect the luminaire to a power supply with an earth conductor (only C-version).
- ▶ Never place additional loads on the luminaire head or the arm.
- ▶ The luminaire must not be covered by a cloth or any similar item while it is in operation.

- ▶ The ventilation openings (where they exist) must be free whenever the luminaire is in operation.
- ▶ The luminaire must not be operated near to external heat sources that exceed the maximum ambient temperature of the luminaire.
- ▶ The luminaire must not be used in environmental conditions other than those for which it is intended.
- ▶ The luminaire must only be used for the intended use described in this document.
- ▶ The manufacturer cannot be held responsible for any injury or damage that is caused as a result of any use other than the intended use or of any failure to comply with safety instructions and warnings.
- ▶ It's recommended to protect the luminaire with an additional earth leakage circuit breaker.

1.4 Warning levels

 DANGER
Indication of hazards that can lead to death or serious injury if measures are disregarded.
 WARNING
Indication of hazards that can lead to injury if measures are disregarded.
CAUTION
Indication of hazards that can lead to damage to property if measures are disregarded.

1.5 Specific mounting instructions

visano 20 C T1


- ▶ **Mounting material is not included.**
- ▶ The ceiling mounting must only be fixed to ceilings that have a concrete stability class B25 (C20/25) or higher.
- ▶ Solid ceiling reinforcement parts must not come into contact with the ceiling mounting. If there is any doubt, an authorized professional must confirm that mounting is possible on the mounting base being used.
- ▶ Borings must be made by specialists in compliance with the boring tolerances for the reinforcement anchor that have been approved by the manufacturer.
- ▶ Screws must be tightened carefully using a torque key that complies with the instructions of the reinforcement anchor manufacturer.

visano 20 P TX

- ▶ **Mounting material is not included.**
- ▶ When using the accessory "wall mount", the installation must be done by a professional.
- ▶ The wall must ensure a firm hold.
- ▶ Use only mounting material which is suitable for the corresponding substrate.

2 SCOPE OF SUPPLY

2.1 visano 20 C T1

- A: 1x ceiling bracket
- B: 1x ceiling cover
- C: 1x end ring
- D: 1x ceiling tube incl. cable
- E: 1x luminaire with adaptor (pre-assembled)

2.2 visano 20 P TX

- ▶ 1x luminaire
- ▶ 1x power cable

3 ASSEMBLY visano 20 C T1

3.1 Load data

Bending moment M_B	25Nm
Vertical weight F_G	90N

3.2 Attach ceiling bracket

⚠ DANGER

Assembly by qualified personnel

- ▶ Assembly must be carried out by qualified personnel only. Lack of appropriate knowledge could be life-threatening.

⚠ DANGER


Life-threatening danger from falling luminaire.

- ▶ The ceiling must be made of solid concrete to guarantee a secure hold.


⚠ DANGER

Electric shocks are life threatening


- ▶ The luminaire must be separated from the mains by an external all-pole switch.


- ▶ Saw ceiling tube to the desired length at the top end using a metal saw and deburr.


- ▶ Remove fastening screw.

- ▶ Insert ceiling bracket 1 into ceiling bracket 2 and drill with $d=9\text{mm}$ 3. Use the existing hole in the ceiling bearing as a guide.

- ▶ Make 4 opposed bore hole markings ($\varnothing 13\text{mm}$), consider the position of the $\varnothing 60\text{mm}$ bore for electrical connections

- ▶ Check clearances

⚠ DANGER

Risk of injury from falling parts

- ▶ The ceiling bracket must be secured to the fastening material, that must be suitable for the corresponding ceiling condition.


⚠ DANGER

Risk of injury from falling parts

- ▶ Mounting must be performed by two people.

CAUTION

Use protective equipment in line with tool manufacturer's instructions


- ▶ Drill the bore holes, then blow them out using bellows

⚠ DANGER

Risk of injury from falling luminaire


- ▶ The ceiling must be made of solid concrete to guarantee a secure hold.
- ▶ Use suitable ceiling anchors.


- ▶ Disassemble the securing element from the ceiling tube (Allen wrench is enclosed, do not lose the threaded pin!)


- ▶ Insert the cable (with loose ends first) and pull up to the stop

- ▶ Pull the luminaire cable of the ceiling tube through the luminaire mount.
- ▶ Insert ceiling tube into luminaire mount
- ▶ Secure using the M8 safety screw and nut 1

- ▶ Hold the ceiling bracket against the ceiling and strike the ceiling anchor with a hammer
- ▶ Tighten fastening in line with manufacturer's data


- ▶ Tighten M8 safety screw and nut 1 (20 Nm / 14.75 lbf ft)
- ▶ Tighten all 4 threaded pins 2 (5 Nm / 3.68 lbf ft)


⚠ DANGER

Mortal danger from electric shock

- ▶ Disconnect or remove fuses before carrying out any work on the power connection

- ▶ Make power connection

- ▶ Connect the secondary cable to the connection point
- ▶ Pay attention to the polarity: **+ = red = brown**
- = blue


- ▶ Push cover and end ring over the ceiling bracket and tighten them (0.5Nm)

3.3 Assemble the luminaire on ceiling tube

⚠ WARNING

The luminaire will fall if fastening is done incorrectly

- ▶ The securing element must be screwed in completely and correctly.

- ▶ Insert the luminaire and the adapter 1 into the ceiling tube 2
- ▶ Insert securing element with threaded pin 3 and screw in to the stop
- ▶ Loosen by one quarter-turn to ensure the function of the turnover axis

4 ASSEMBLY visano 20 P TX

4.1 Load data

Bending moment M_B	25Nm
Vertical weight F_G	20N

4.2 Assembly

- ▶ The luminaire are equipped with an adapter pin. The luminaire must be positioned in one of the accessories described in chapter 9.

5 OPERATION

5.1 visano 20 C T1

⚠ DANGER

Electric shocks are life-threatening

- ▶ Only connect to a mains power supply with a protective conductor

5.2 visano 20-2

⚠ DANGER


Electric shocks are life-threatening

- ▶ Do not insert any power cable that is damaged.
- ▶ If there is any sign of damage to the power cable, replace it immediately with a new one.
- ▶ The connection voltage and frequency must match the data on the type plate.


⚠ WARNING

Risk of eye damage


- ▶ Never look directly into the beam of light.

- ▶ Insert power cable to the luminaire
- ▶ Connect the cable with the mains
- ▶ Perform a function test before every use: All LED's in the light cone must illuminate.


- ▶ Switch the LED module on/off


- ▶ Set the desired colour temperature according to the following table:

COLOUR	COL. TEMP.	RECOMMENDATION
blue	cold white	Ideal for identifying individual vessels for wound follow-up treatment and small procedures
orange	warm white	Ideal for various types of skin examinations

- ▶ Set the desired Illuminance in 25% dimmer steps

5.2 Transport position

- ▶ Transport the luminaire only in this position.

6 DISINFECTION CLEANING

DANGER

Electric shocks are life-threatening

- ▶ Remove the plug from the mains before the disinfection cleaning.

CAUTION

Damage to property caused by incorrect cleaning

- ▶ Use only disinfectants that will not affect the functional capability of the luminaire.
- ▶ No cleaning agents containing alcohol, solvents, chlorine or abrasive products must be used. Those agents can cause damages in the plastic parts.
- ▶ The cleaning and disinfection agents must be suitable and approved for plastics like PC, PMMA, PA and ABS.
- ▶ Concentrated disinfectants may damage the cover.
- ▶ Pay attention to the specifications from the data sheet of the agent for concentration and residence time.
- ▶ Unsuitable cloths can cause scratches.

NOTE

Dirt decreases luminosity

- ▶ Clean regularly to keep the screen clean.
- ▶ Disinfection only by wiping is permitted.

- ▶ Clean the PA screen with a not abrasive cleaning cloth and a suitable cleaning agent.

CAUTION

To minimize the risk of disease transmission, in addition to complying with this instruction for use, you must also comply with the applicable occupational health and safety regulations and the requirements of national bodies with responsibilities for hygiene and disinfection.

7 SAFETY INSPECTIONS

DANGER

Electric shocks are life-threatening

- ▶ Remove the plug from the mains.
- ▶ The connecting cable must be checked for damage at least once per year.

NOTE

- ▶ Maintenance and repairs must only be carried out by qualified electricians.
- ▶ The corresponding user profile is described in Chapter 1, Safety instructions

8 DEMOUNTING

DANGER

Electric shocks are life-threatening

- ▶ Remove the plug from the mains before demounting.

8.1 Disposal

Do not dispose the luminaire in household refuse. Dispose the luminaire at a disposal point in accordance with local regulations or take them to a dealer that provides an appropriate disposal service. Cut off the cable directly on the housing.

The products listed above are more than 95% recyclable. The luminaires have been constructed to be compatible with recycling so that a high proportion of the materials used in these products can be recycled or converted into energy after the end of service life. They contain no materials that are dangerous or that need to be monitored.

9 ACCESSORIES

WARNING

The luminaire will fall if fastening is done incorrectly

- ▶ Use the extension arm (D13.363.000) only in combination with the wall mount (new) (D15.435.000).

Extension arm (Art. Nr.: D13.363.000 - pure white)

Table clamp (Art. Nr.: D14.228.000 - black)

Universal mount (Art. Nr.: D13.430.000 - pure white)

Rail fastening (Art. Nr.: D13.269.000)

Wall mount (Art. Nr.: D13.231.000 - pure white)

Wall mount (new) (Art. Nr.: D15.435.000 - pure white)

Roller stand RL (Art. Nr.: D14.022.000)

10 ADDITIONAL INSTRUCTIONS

The luminaire itself is maintenance-free.

You can obtain additional documents about this product from the manufacturer upon request.

These luminaires do not affect any other devices.

To save energy, the luminaire should only be switched on when actually needed.

11 TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profiles
The luminaire does not light	Contact problem	Switch on again	All
The luminaire does not light	LED module defective	Contact manufacturer service	By manufacturer service only
The luminaire does not light	No mains voltage	Check mains voltage, check all connections	Electrician

12 TECHNICAL DATA

Electrical data:	
Rated input voltage	100 - 240VAC
Frequency range	50/60 Hz
Power consumption	12-15W (14-20VA)
Input current	0.08-0.14A
Power factor	0.75-0.86
Power supply secondary side	24VDC
Photometric values:	
Central illuminance E_c at 0.5m (1.64 feet) distance	60'000 lx *
Light field diameter d_{10} at 0.5m (1.64 feet) distance	$\varnothing = 21\text{ cm (8.3")}$ *
Light field diameter d_{50} at 0.5m (1.64 feet) distance	$\varnothing = 12\text{ cm (4.7")}$ *
Colour temperature	4500K / 3500K *
Colour rendering index R_a	>95
Colour rendering index R_9	>90
Total irradiance E_e at maximum intensity	<230 W/m ²
Ratio irradiance E_e to illuminance E_c	3.83 mW/m ² /lx
	* -10% / +20% tolerance
Ambient conditions for transport, storage and operation:	
Ambient temperature (storage and transport)	-20°C to +70°C (-4°F to +158°F)
Ambient temperature (operation)	+10°C to +35°C (+50°F to +95°F)
Relative humidity (non-condensing) (storage and transport)	max. 90%
Relative humidity (non-condensing) (operation)	max. 75%
Weight:	
Luminaire head	0.6kg (1.3 lb)
visano 20 C T1	9kg (19.8 lb)
visano 20 P TX	2kg (4.4 lb)
Mode of operation:	
Mode of operation	Continuous operation

Classification:	
visano 20	Protection class I (C-version) Protection class II (P-version)
Degree of protection according to IEC 60529	IP 20
Luminaire head (only in horizontal position)	IP 43
Classification according to directive 93/42/EWG – annex IX (medical device class)	Class I
U.S. FDA Device Class	Class I
Electrical safety test and EMC according to:	IEC 60601-1 IEC 60601-2-41 IEC 60601-1-2
GMDN-Code	12276
Service life of the light source:	
Service life LED	30'000h (L70/B50)

13 ELECTROMAGNETIC COMPATIBILITY (EMC)

Electromagnetic emission guidelines		
The medical device is intended for operation in an electromagnetic environment such as the one specified below. The user must ensure that it is operated in such an environment.		
Emissions	In accordance	Electromagnetic environment
High-frequency emissions (CISPR 11)	Group 1	The medical device uses HF energy exclusively for its internal function. This means that its HF emissions are very low, and it is unlikely that adjacent electronic devices will receive interference.
High-frequency emissions (CISPR 11)	Class B	The medical device is intended for use in all facilities, including residential buildings and facilities that are directly connected (without a transformer) to the same low voltage network as residential buildings.
Emissions from harmonics (IEC 61000-3-2)	Class A	
Emissions from voltage fluctuations/flicker (IEC 61000-3-3)	Compliance	

Electromagnetic interference resistance guidelines			
The medical device is intended for operation in an electromagnetic environment such as the one specified below. The user must ensure that it is operated in such an environment.			
Interference resistance to	IEC 60601-1-2 testing level	Conformity level of the medical device	Electromagnetic environment
Electrostatic discharge (ESD) (IEC 61000-4-2)	Contact discharge: ± 6 kV Air discharge: ± 8 kV	± 6 kV ± 8 kV	Floors made of timber, concrete or ceramics tiles are preferred. Where synthetic floor covering is used the relative humidity should be at least 30%.
Fast transients/ electrical disturbance variables/ bursts (IEC 61000-4-4)	Power supply cables: ± 2 kV Longer input and output cables: ± 1 kV	± 2 kV Not applicable	The quality of the mains power supply should match that of a typical business or hospital environment.
Surge voltages/surges (IEC 61000-4-5)	±1 kV phase-to-phase voltage ±2 kV phase-to-earth voltage	±2 kV ±1 kV	The quality of the mains power supply should match that of a typical business or hospital environment.
Magnetic field in the power supply frequency (50/60 Hz) (IEC 61000-4-8)	3 A/m	3 A/m	In proximity of this medical device, do not operate devices with unusually strong network-frequency magnetic fields (transformer stations etc.) .
Voltage dips and short interruptions to the power voltage (IEC 61000-4-11)	Cut >95 %, 0.5 periods Cut 60 %, 5 periods Cut 30 %, 25 periods Cut >95 %, 5 seconds	Cut >95 %, 0.5 periods Cut 60 %, 5 periods Cut 30 %, 25 periods Cut >95 %, 5 seconds	The supply voltage quality should be that of a typical business or hospital environment. If the user requires continued function during any interruption of the energy supply system, we recommend that the medical device be powered by an uninterrupted power supply or a battery.
Radiated HF disturbance variable (IEC 61000-4-3)	3 V/m 80 MHz - 2.5 GHz	3 V/m	Recommended minimum distance of portable and mobile radio devices of PEIRP radiated power to the medical device, including its cables: $d = 1.17\sqrt{P}$
Conducted HF disturbance values (IEC 61000-4-6)	$3 V_{\text{rms-value}}$ 150 kHz - 80 MHz:	$3 V_{\text{rms-value}}$	Recommended minimum distance of portable and mobile radio devices of PEIRP radiated power to the medical device, including its cables: 80MHz - 800MHz: $d = 1.17\sqrt{P}$ 800MHz - 2.5GHz: $d = 2.33\sqrt{P}$
d = recommended safe distance [m], P = Rated Power of the sender [W]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range.			
Interference may occur in the vicinity of equipment marked with: 			

Recommended safe distances from portable and mobile HF communications equipment		
Rated Power of the sender [W]	150 kHz - 800 MHz	800 MHz - 2.5 GHz
0.01	0.12 m (0.39 ft)	0.23 m (0.75 ft)
0.1	0.37 m (1.21 ft)	0.74 m (2.43 ft)
1	1.17 m (3.84 ft)	2.33 m (7.64 ft)
10	3.7 m (12.14 ft)	7.37 m (24.18 ft)
50	8.27 m (27.13 ft)	16.48 m (54.07 ft)
100	11.7 m (38.39 ft)	23.3 m (76.44 ft)

OWN DISTRIBUTORS

GERMANY

Herbert Waldmann GmbH & Co. KG
Postfach 5062
78057 Villingen-Schwenningen
Germany
Telephone +49 7720 601 0
Telephone +49 7720 601 100 (Sales)
Fax +49 7720 601 290
www.waldmann.com
sales.germany@waldmann.com

SWITZERLAND

Waldmann Lichttechnik GmbH
Benkenstrasse 57
5024 Küttigen
Switzerland
Telephone +41 62 839 12 12
Fax +41 62 839 12 99
www.waldmann.com
info-ch@waldmann.com

AUSTRIA

Waldmann Lichttechnik Ges.m.b.H
Gewerbepark Wagram 7
4061 Pasching/Linz
Austria
Telephone +43 7229 67 400
Fax +43 7229 67 444
www.waldmann.com
info-at@waldmann.com

FRANCE

Waldmann Eclairage S.A.S
Z.I. - Rue de l'Embranchement
67116 Reichstett
France
Telephone +33 3 8820 95 88
Fax +33 3 8820 95 68
www.waldmann.com
info-fr@waldmann.com

ITALY

Waldmann Illuminotecnica S.r.l.
Via della Pace, 18 A
20098 San Giuliano Milanese (MI)
Italy
Telephone +39 02 98 24 90 24
Fax +39 02 98 24 63 78
www.waldmann.com
info-it@waldmann.com

NETHERLANDS

Waldmann BV
Lingewei 19
4004 LK Tiel
Netherlands
Telephone +31 344 631 019
Fax +31 344 627 856
www.waldmann.com
info-nl@waldmann.com

SWEDEN

Waldmann Ljusteknik AB
Skebokvarnsvägen 370
124 50 Bandhagen
Sweden
Telephone +46 8 990 350
Fax +46 8 991 609
www.waldmann.com
info-se@waldmann.com

USA

Waldmann Lighting Company
9, W. Century Drive
Wheeling, Illinois 60090
USA
Telephone +1 847 520 1060
Fax +1 847 520 1730
www.waldmannlighting.com
waldmann@waldmannlighting.com

SINGAPORE

Waldmann Lighting Singapore Pte. Ltd.
22 Cross Street, #02-50/51
South Bridge Court
Singapore 048421
Singapore
Telephone +65 6275 8300
Fax +65 6275 8377
www.waldmann.com
sales-sg@waldmann.com

CHINA

Waldmann Lighting (Shanghai) Co., Ltd.
Part A11a, No. Five Normative Workshop
199 Changjian Road, Baoshan
Shanghai, P.R.C. 200949
China
Telephone +86 21 5169 1799
Fax +86 21 3385 0032
www.waldmann.com
info@waldmann.com.cn

INDIA

Waldmann Lighting PVT. Ltd.
Plot No. 52
Udyog Vihar
Phase-VI, Sector-37
GURGAON-122011, Haryana
INDIA
Telephone +91 124 412 1600
Fax +91 124 412 1611
www.waldmann.com
sales-in@waldmann.com

Diese Montage- und Gebrauchsanweisung dient ausschliesslich der Kundeninformation und wird nur auf Kundenanforderung aktualisiert oder ausgetauscht.
These installation and operating instructions are for customer information only and will only be updated or replaced upon request by the customer.
Ces instructions d'installation et de fonctionnement sont destinées au client uniquement et ne seront mises à jour ou remplacées uniquement sur demande du client.
Le presenti istruzioni per il montaggio e l'uso servono esclusivamente come informazione per il cliente e sono aggiornate o sostituite solo su richiesta del cliente.
Deze montage- en gebruiksaanwijzing is uitsluitend bestemd als informatie voor de klant en wordt alleen op verzoek van de klant bijgewerkt of vervangen.
Este manual de instrucciones de montaje y utilización sirve exclusivamente para informar al cliente y solo se actualiza o sustituye a petición del cliente.
これらの取り付けと操作の手引きはお客様向けに作成したものです。アップデートや交換は要求してください。

Änderungen vorbehalten • gedruckt in der Schweiz | Subject to changes • printed in Switzerland | Peut être modifié • imprimé en Suisse
Con riserva di modifiche • Stampato in Svizzera | Wijzigingen voorbehouden • Gedrukt in Duitsland | Salvo modificaciones • impreso en Alemania

Derungs Licht AG
Hofmattstrasse 12
CH-9200 Gossau SG
Switzerland
Telephone +41 71 388 11 66
Fax +41 71 388 11 77

Further distribution partners you find at: www.derungslicht.com