



TRIANGO 60-1 F



Fitted with	LED - Energy efficiency category A+
Illuminance field 100 cm	ø 18 cm
Illuminance 100 cm	60 000 lx
Colour rendering index (Ra)	> 94
Work equipment	Electronic converter
Voltage	100-240 V; 50/60 Hz
Mains lead	Direct power supply
Luminaire body	Polyamid (PA)
Weight (net)	19 kg
Power consumption	< 30 W
Usage	Keypad at cardanic
Balance of articulation joints	Spring
Class of protection	I
Accessory	Sterilizable handle, handle cover
Norms	EN 60601-1, EN 60601-2-41
Special features	Stepless dimmable, diamond optics, sterilizable handle
Fastening	Roller stand
System of protection	IP 43 (luminaire head)
Medical device class	I according to MDR

	Lighting field / 1,0 m fixed value 18cm	lighting field / 1,0 m changeable value 18 / 23 / 25cm	Colour temperature fixed value 4700K	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	Order No.
TRIANGO 60-1 F CH	X		X		X	D16.262.000-00813744
TRIANGO 60-1 F Schuko	X		X		X	D16.259.000-00813743
TRIANGO 60-1 F UK	X		X		X	D16.271.000-00813747
TRIANGO 60-1 F AUS	X		X		X	D16.268.000-00813746
TRIANGO 60-1 F US	X		X		X	D16.265.000-00813745

Light distribution



Technical drawing

