


Treatment luminaire

TRIANGO Endo 100-1 C D15972000-00803704



fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm
illuminance @ 100 cm	100 000 lx
color rendering index (CRI)	> 95
colour temperatur	4300 K
work equipment	electronic converter
mains lead	direct power supply
luminaire body	polyamide (PA)
weight (net)	15.9 kg
Power consumption	40 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle, adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	4 dimming levels, Diamond optics, Endo Mode, High-quality arm system, sterilisable grip, maintenance free
Fastening	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
System of protection	IP 20 (43)
order number	D15972000 - pure white

light distribution

