

Perimetry Trial Lens Set-Plus Cylinder Set 68 PCS
(leatheroid case, silver metal rim)

SPHERE LENS								CYLINDER LENS	
CONVEX(+)				CONCAVE(-)				+Cylinder	
list	pcs	list	pcs	list	pcs	list	pcs	list	pcs
0.25	1	4.50	1	0.25	1	4.50	1	0.25	1
0.50	1	5.00	1	0.50	1	5.00	1	0.50	1
0.75	1	5.50	1	0.75	1	5.50	1	0.75	1
1.00	1	6.00	1	1.00	1	6.00	1	1.00	1
1.25	1	6.50	1	1.25	1	6.50	1	1.50	1
1.50	1	7.00	1	1.50	1	7.00	1	2.00	1
1.75	1	8.00	1	1.75	1	8.00	1	2.50	1
2.00	1	9.00	1	2.00	1	9.00	1	3.00	1
2.25	1	10.00	1	2.25	1	10.00	1	3.50	1
2.50	1	11.00	1	2.50	1	11.00	1	4.00	1
2.75	1	12.00	1	2.75	1	12.00	1	4.50	1
3.00	1	13.00	1	3.00	1	13.00	1	5.00	1
3.50	1	14.00	1	3.50	1	14.00	1		
4.00	1	16.00	1	4.00	1	16.00	1		

Perimetry Trial Lens Set-Minus Cylinder Set 68 PCS
(leatheroid case, silver metal rim)

SPHERE LENS								CYLINDER LENS	
CONVEX(+)				CONCAVE(-)				-Cylinder	
list	pcs	list	pcs	list	pcs	list	pcs	list	pcs
0.25	1	4.50	1	0.25	1	4.50	1	0.25	1
0.50	1	5.00	1	0.50	1	5.00	1	0.50	1
0.75	1	5.50	1	0.75	1	5.50	1	0.75	1
1.00	1	6.00	1	1.00	1	6.00	1	1.00	1
1.25	1	6.50	1	1.25	1	6.50	1	1.50	1
1.50	1	7.00	1	1.50	1	7.00	1	2.00	1
1.75	1	8.00	1	1.75	1	8.00	1	2.50	1
2.00	1	9.00	1	2.00	1	9.00	1	3.00	1
2.25	1	10.00	1	2.25	1	10.00	1	3.50	1
2.50	1	11.00	1	2.50	1	11.00	1	4.00	1
2.75	1	12.00	1	2.75	1	12.00	1	4.50	1
3.00	1	13.00	1	3.00	1	13.00	1	5.00	1
3.50	1	14.00	1	3.50	1	14.00	1		
4.00	1	16.00	1	4.00	1	16.00	1		