

UBM Plus®

Item Number: 24-6300

Accutome revolutionized the world of ultrasound imaging with the launch of the B-Scan Plus®. We proudly introduce the new portable, high definition UBM Plus. The state-of-the-art probe plugs into a laptop or PC. The anterior segment images are recorded in a 34 second film loop and easily transfer to EMR. Upgradeable software protects your investment.

Features:

- High definition imaging - Detailed structure definition including cornea, iris, ciliary body, zonules, crystalline lens and intraocular lens as well as pathologies.
- State-of-the-art probe design - Sharper, more focused images due to the elimination of signal loss.
- Unsurpassed data analysis - Contains measuring tools. Useful for measuring sulcus-to-sulcus, anterior chamber depth, positioning of intraocular lenses and filtration angle of the eye.
- Portable - Plug probe into any laptop or PC and go.
- Never miss an image - Unlimited 34 second film loops.
- Share information easily - Adaptable document transfer via EMR, email or printer.
- Protect your investment - Easily upgrade software.
- Master in minutes - User-friendly interface makes user manuals and quick start guides unnecessary.

