



TRIANGO 60 C



Fitted with	LED - Energy efficiency category A+
Illuminance field 100 cm	ø 18 cm
Illuminance 100 cm	60 000 lx
Colour rendering index (Ra)	> 94
Work equipment	Electronic converter
Voltage	100-240 V; 50/60 Hz
Mains lead	Direct power supply
Luminaire body	Polyamid (PA)
Weight (net)	13 kg
Power consumption	< 30 W
Usage	Keypad at cardanic
Balance of articulation joints	Spring
Class of protection	I
Accessory	Sterilizable handle, handle cover, adapter for suspended ceiling
Norms	EN 60601-1, EN 60601-2-41
Special features	Stepless dimmable, diamond optics, sterilizable handle
Fastening	Ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
System of protection	IP 43 (luminaire head)
Medical device class	I according to MDR

	lighting field / 1,0 m fixed value 18cm	lighting field / 1,0 m changeable value 18 / 23 / 25cm	Colour temperature fixed value 4700K	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	Order No.
TRIANGO 60-1 C	X		X		X	D16.250.000-00813741
TRIANGO 60-1 C US	X		X		X	D16.253.000-00814312
TRIANGO 60-3 C	X			X	X	D16.251.000-00814310
TRIANGO 60-3 C US	X			X	X	D16.254.000-00814313
TRIANGO Fokus 60-1 C		X	X		X	D16.252.000-00814311
TRIANGO Fokus 60-1 C US		X	X		X	D16.255.000-00814314

Light distribution



Technical drawing

