

Treatment luminaire

TRIANGO 80-1 C D16131000-00808951



Illuminancefield @ 100 cm	Ø 18 cm
illuminance @ 100 cm	Ø 80 000 lx (with dimming 32 000 lx to 80 000 lx)
color rendering index (CRI)	> 95
colour temperatur	4500 K
Depth of illumination L1+L2	1.24 m
work equipment	electronic converter
mains lead	direct power supply
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	14.0 kg
Power consumption	30 W
usage	rotary switch on the luminaire head single-pole secondary, continuously dimmable
balance of articulation joints	spring (weight compensation)
class of protection	I
Accessory	sterilizable handle, adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	sterilisable grip, continuously dimmable, maintenance free
Fastening	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm
System of protection	IP 20 (43)
order number	D16131000 - pure white

light distribution

