

Minor surgical light

TRIANGO Fokus 100-3 C DUO D16202000-00811330

Derungs



MEDICAL LIGHTING



fitted with	LED - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 18 cm (each)
illuminance @ 100 cm	100 000 lx (each)
color rendering index (CRI)	> 95 (each)
colour temperatur	3700/4300/4700 K (each)
work equipment	electronic converter
mains lead	direct power supply (2x)
luminaire body	polyamide (PA)
weight (net)	30.0 kg
Power consumption	80 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	3 color shades, 4 dimming levels, Diamond optics, Focusable, High-quality arm system, sterilisable grip, maintenance free
Fastening	ceiling mount, when mounting ceiling tube could be cut to length 25-120 cm, intermediate ceiling adapter (on request)
System of protection	IP 20 (43)
order number	D16202000 - pure white

light distribution

