



TRIANGO 80-1 F Battery



Fitted with	LED - Energy efficiency category A+ ø 18 cm
Illuminance field 100 cm	80 000 lx
Illuminance 100 cm	> 95
Colour rendering index (Ra)	Electronic converter
Work equipment	100-240 V; 50/60 Hz
Voltage	Direct power supply
Mains lead	Polyamid (PA)
Luminaire body	32 kg
Weight (net)	30 W
Power consumption	Rotary switch on the luminaire head
Usage	Spring
Balance of articulation joints	I
Class of protection	Sterilizable handle, handle cover
Accessory	EN 60601-1, EN 60601-2-41
Norms	Continuously dimmable, diamond optics, sterilizable handle
Special features	Roller stand
Fastening	IP 43 (luminaire head)
System of protection	IP 20 (UPSmed pack)
Medical device class	I according to MDR
Battery cycle life up to 70%	1.200 cycles
Battery type	NP17-12I Industrial VRLA Battery

	Lighting field / 1,0 m fixed value 18 cm	lighting field / 1,0 m changeable value 18 - 25 cm	Colour temperature fixed value 4500K	Colour temperature changeable value 3700K / 4300K / 4700K	Dimmable	UPSmed Pack	Order No.
TRIANGO 80-1 F Battery Schuko	X		X		X	X	D16.209.000
TRIANGO 80-1 F Battery CH	X		X		X	X	D16.218.000
TRIANGO 80-1 F Battery UK	X		X		X	X	D16.219.000
TRIANGO 80-1 F Battery US	X		X		X	X	D16.220.000
TRIANGO 80-1 F Battery AUS	X		X		X	X	D16.221.000

Light distribution



Technical drawing

